


MEASUREMENT/FITTING GUIDE

Patient Name: _____ Date: _____

FITTING TIPS

- If the calf measurement is one size larger than the ankle on the fitting chart, the size should be based on the measurement of the calf.
- If the ankle measurement is one size larger than the calf on the fitting chart, the size should be based on the measurement of the ankle.
- If the difference between the calf and the ankle is more than one size in either direction, the garment should not be used and a custom prescribed compression wear should be utilized.
- If both the ankle and calf measurements fall at the mid-point between sizes, select the smaller size to avoid quickly sizing out of the garment. As fluid is dispersed, the smaller size will allow for longer use.



TIBIAL TUBEROSITY - Bony "Bump" below knee

CALF MEASUREMENT - Taken at widest part of the calf muscle

LENGTH MEASUREMENT - Taken from floor to Tibial Tuberosity

ANKLE MEASUREMENT - Taken just above the ankle

EXTREMIT-EASE® COMPRESSION GARMENT

LEG CIRCUMFERENCE MEASUREMENTS

LENGTH MEASUREMENTS

	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE	REG	TALL
A = ANKLE	21 - 25 cm 8.5 - 10 in	25 - 30 cm 10 - 11.5 in	30 - 33 cm 11.5 - 13 in	33 - 37 cm 13 - 14.5 in	37 - 41 cm 14.5 - 16 in	37 - 44 cm 14.5 - 17.5 in	44 - 51 cm 17.5 - 20 in
C = CALF	37 - 42 cm 14.5 - 16.5 in	42 - 47 cm 16.5 - 18.5 in	47 - 52 cm 18.5 - 20.5 in	52 - 57 cm 20.5 - 22.5 in	57 - 62 cm 22.5 - 24.5 in		

EXTREMIT-EASE® COMPRESSION GARMENT/LINER ORDER:

One Garment Liner and mesh laundry bag are included with each EXTREMIT-EASE® COMPRESSION GARMENT

Indicate number/size needed:

NO product substitutions

REGULAR

- ___ XS (E10140)
- ___ S (E10141)
- ___ M (E10142)
- ___ L (E10143)
- ___ XL (E10144)

TALL

- ___ XS (E20140)
- ___ S (E20141)
- ___ M (E20142)
- ___ L (E20143)
- ___ XL (E20144)

ADDITIONAL GARMENT LINERS

- ___ QTY - M (E90002)
- ___ QTY - L (E90003)

EXTREMIT-EASE® LINER

	MEDIUM	LARGE
MALE SHOE SIZE	7 - 9	10 - 12
FEMALE SHOE SIZE	8 - 10	11 - 13

SPECIFIC INSTRUCTIONS: _____