



PARTNERS IN THE SCIENCE &
MANAGEMENT OF WOUND CARE



An abstract background featuring a blue and white geometric pattern of interconnected triangles and lines, creating a sense of depth and movement. The pattern is more dense and vibrant in the upper left and center, fading into a lighter blue and white gradient towards the bottom right.

CONTACT US WITH ANY QUESTIONS!

AMERX Health Care Corporation

(800) 448-9599 • info@amexhc.com • AMERXHC.com

TABLE OF CONTENTS

TURN-KEY DME® Wound Care Program2

AMERX® Wound Care Kits4

EXTREMIT-EASE® Compression Garments & Liners6

AMERIGEL® Advanced Skin & Wound Care8

HELIX3® Bioactive Collagen14

AMERX® Dressings17

About AMERX Health Care

Founded in 1993 by John C. Anderson, AMERX Health Care Corporation set out to develop quality, trusted and economical wound care products for at-risk patients. Our flagship product line, AMERIGEL®, formed the foundation for what would become the standard of excellence at AMERX. Top rankings and a reputation for positive results encouraged AMERX to develop a comprehensive suite of wound care products including HELIX3® Bioactive Collagen, AMERX® Wound Care Kits & Dressings, EXTREMIT-EASE® Compression Garment and our newest addition, ADVANTAGEN® Surgical Collagen.

Our high-quality, affordable products address today’s wound and skin care challenges – providing health care facilities and private practice physicians the tools for success while helping patients live a healthy life.



TURN-KEY DME[®] WOUND CARE PROGRAM

Don't leave your DME Wound Care Program to chance!

AMERX Health Care's Turn-Key DME[®] Program provides over 25 years of PROVEN success in delivering the support, training, and tools necessary to drive a successful DME Wound Care Program.

Experience a True Partnership with Ongoing Support

Your success is our first priority! AMERX support starts with our first-in-class onboarding program. Your dedicated AMERX team member will connect you to top coding and billing consultants and personalize our Turn-Key DME program to suit your practice needs. Once completed, your AMERX representative will work directly with your office manager, Medical Assistants and billing team to provide ongoing support and ensure your program is exceeding your expectations.

Enhance Your Business Revenues with Billing and Coding Consultations from Nationally-Known Experts

AMERX support extends beyond our knowledgeable team members to include upfront and ongoing assistance from nationally recognized experts within the fields of Coding and Billing. This team of experts continuously monitors for changes/updates to Medicare policies throughout the year to provide your practice with up-to-date information, so you can focus on your patient's needs.

Discover the Value of Turn-Key DME! (800) 448-9599 | join@amerxhc.com



Enjoy Access to a Full Suite of High-Quality Wound Care Products

Your options are never limited! Purchase only the products you want to dispense direct from your office and order any other products from AMERX Health Care's expansive line of wound care solutions to be delivered direct to the patient.

PATIENT DIRECT

Zero-Inventory Dispensing Program

- ▶ Drive results with access to the full suite of AMERX wound care products.
- ▶ Product delivered direct to patient's doorstep.
- ▶ Verified insurance coverage and accurate claims submission.
- ▶ NO upfront costs or minimums.

Experience the convenience of access to a full suite of AMERX products delivered direct to your patients.

DOCTOR DIRECT

In-Office Dispensing Program

- ▶ Provide the convenience today's patients expect.
- ▶ Patients leave office with product in-hand.
- ▶ Improve patient compliance with brands you trust.
- ▶ Volume discount pricing.

Experience the cost savings of dispensing DME wound care supplies direct to your patients.

Top Reasons to Dispense AMERX Products

BOOST POSITIVE OUTCOMES - Trusted, FDA-cleared, PDAC-approved products with proven results for more than 25 years.

INCREASE PATIENT COMPLIANCE - Studies show 60-70% higher patient compliance rates when products are supplied direct to the patient.

ENHANCE PATIENT SATISFACTION - Save your patients time and money by helping them avoid additional trips to the pharmacy.

WOUND CARE KITS WITH BORDERED GAUZE

INCREASE Patient Wound Care Compliance

Everything needed for your patients to perform proper dressing changes at home including:

- ✓ Primary Dressing
- ✓ Secondary Dressing
- ✓ Saline Wound Wash
- ✓ Supportive Dressing

Available in a variety of primary and secondary dressing options based on individual wound needs.

Reimbursable through Medicare and most private insurance plans.



Primary Dressing	AMERX Bordered Gauze		Catalog / SKU #	HCPCS Codes
	Border Size	Pad Size		
30-Day Kits				
HELIX3-CP® Collagen Powder (1g.)	2x2	1x1	160419	A6010, A6219, A6216
HELIX3-CP® Collagen Powder (1g.)	4x4	2x2	160416	A6010, A6219, A6216
HELIX3-CM® Collagen Matrix (2x2)	4x4	2x2	160410	A6021, A6219, A6216
AMERX® Calcium Alginate (2x2)	4x4	2x2	160620	A6196, A6219, A6216
AMERX® Foam Dressing (2x2)	4x4	2x2	160517	A6209, A6219, A6216
AMERIGEL® Hydrogel Wound Dressing (3oz.)	4x4	2x2	160204	A6248, A6219, A6216
15-Day Kits				
HELIX3-CP® Collagen Powder (1g.)	2x2	1x1	160418	A6010, A6219, A6216
HELIX3-CP® Collagen Powder (1g.)	4x4	2x2	160415	A6010, A6219, A6216
HELIX3-CM® Collagen Matrix (2x2)	4x4	2x2	160409	A6021, A6219, A6216
AMERX® Calcium Alginate (2x2)	4x4	2x2	160619	A6196, A6219, A6216
5-Day Kits				
HELIX3-CP® Collagen Powder (1g.)	2x2	1x1	160417	A6010, A6219, A6216



NEW!

WOUND CARE KITS WITH ROLLED GAUZE



Primary Dressing	Secondary Dressings			Catalog / SKU #	HCPCS Codes
	Rolled Gauze	Sterile Gauze Pads	Paper Tape		
30-Day Kits					
HELIX3-CP® Collagen Powder (1g.)	3 in.	4x4	1 in.	160431	A6010, A6446, A6402, A4450
HELIX3-CM® Collagen Matrix (2x2)	3 in.	4x4	1 in.	160433	A6021, A6446, A6402, A4450
AMERX® Calcium Alginate (2x2)	3 in.	4x4	1 in.	160635	A6196, A6446, A6402, A4450
AMERX® Foam Dressing (2x2)	3 in.	4x4	1 in.	160536	A6209, A6446, A6402, A4450
AMERIGEL® Hydrogel Wound Dressing (3oz.)	3 in.	4x4	1 in.	160237	A6248, A6446, A6402, A4450
15-Day Kits					
HELIX3-CP® Collagen Powder (1g.)	3 in.	4x4	1 in.	160430	A6010, A6446, A6402, A4450
HELIX3-CM® Collagen Matrix (2x2)	3 in.	4x4	1 in.	160432	A6021, A6446, A6402, A4450
AMERX® Calcium Alginate (2x2)	3 in.	4x4	1 in.	160634	A6196, A6446, A6402, A4450

Contact AMERX for HCPCS Coding Guidance.

COMPRESSION THERAPY



Provides 30-50 mmHG Compression

INDICATIONS:

Management of mild to moderate swelling from lymphedema, chronic venous insufficiency, acute post-operative edema and other swelling disorders.

Gradient Compression Wrap

No Natural Rubber Latex

Includes Garment, Liner, and Laundry Bag

Reimbursable through Medicare and most private insurance plans.



EXTREMIT-EASE®
COMPRESSION GARMENT

- Easy for patients to independently apply and adjust.
- Made with lightweight, air permeable fabric for hours of comfortable wear.
- Reduces pistoning, common with layered compression.
- Machine washable with provided mesh laundry bag.

Color	Length	Size	Catalog / SKU#	HCP Code
Black	Regular	XS – XXL	E10140 – E10145	A6545
Black	Tall	XS – XXL	E20140 – E20145	A6545
Tan	Regular	XS – XXL	E10240 – E10245	A6545
Tan	Tall	XS – XXL	E20240 – E20245	A6545

EXTREMIT-EASE® GARMENT LINER

- Combines the benefit of mild mid-foot to ankle compression stocking with the comfort of a low-profile liner.
- Worn independently or in conjunction with EXTREMIT-EASE Compression Garments.
- Low-profile foot compression that can be worn with dress shoes.



**Provides 10-20 mmHg
Foot Compression**

RECOMMENDED FOR:
Mild to moderate edema

Easy to Apply
Machine Washable
No Natural Rubber Latex

Color	Size	Catalog / SKU#	HCPCS Code
Gray	S – XXL	E90011 – E90015	N/A
Tan	S – XXL	E90001 – E90005	N/A

Podiatrist Recommended for Diabetics' Skin Relief

RECOMMENDED FOR:

Rehydrating and soothing dry, irritated skin, softening and exfoliating corns and calluses.

Moisturizing and protecting dry skin associated with diabetes.

Minimizing the appearance of bruising (hemosiderin staining).

Venous stasis dermatitis.

Fragrance, Dye and Paraben-FREE



**100% Direct-to-Patient
Money Back Guarantee**



AMERIGEL[®] CARE LOTION

- Advanced micro-particle technology targets conditions beneath the skin's surface.
- Fast absorbing formula is non-greasy, non-comedogenic and hypoallergenic.
- Safe to use between the toes without causing maceration.
- Ideal for a regular skin routine with compression therapy.

Size	Catalog / SKU#	HCPCS Code
6 oz. Tube	A30106	N/A
16 oz. Pump Bottle	A30116	N/A

AMERIGEL® BARRIER LOTION

- Invisible shield for advanced skin protection.
- Contains 5% Dimethicone for premium skin protection.
- Dual action skin protectant shields and lubricates skin.
- Effortless application and removal.
- Ideal for use as skin protectant under compression garments and prostheses.



**Physician Recommended for
Use with Compression Garments**

INDICATIONS:

Helps prevent and temporarily protect chafed, chapped or cracked skin. Helps prevent and protect from the drying effects of wind and cold weather. Helps protect from wetness.

Non-greasy

Dries clear



**100% Direct-to-Patient
Money Back Guarantee**

Size	Catalog / SKU#	HCPSC Code
6 oz. Tube	A30206	N/A
16 oz. Pump Bottle	A30216	N/A

Available Exclusively
Through Physician Offices

RECOMMENDED FOR
USE FOLLOWING:

Nail surgery, skin biopsy,
wart removal, wound
debridement, laser surgery
and surgical procedures.

KIT CONTAINS:

- 1 oz. AMERIGEL
Hydrogel Dressing
- 4 oz. AMERIGEL
Saline Wash
- 30 2x2 Gauze Pads

3 Secondary Dressings
To Choose From!

- Fabric Bandage
- Knuckle Bandage
- Flex Tape



AMERIGEL®
POST-OP
SURGICAL KIT

- More cost effective and convenient for your patients, saving time and money while avoiding trips to the pharmacy.
- Patient compliance rates are **60-70%** higher when patients leave the practice with product in-hand vs. filling a prescription.
- Contains the #1 Rated Topical Dressing for Wound/Ulcer Treatment and Treatment following Nail Surgery since 2003.*

*Podiatry Management Annual/Physicians Surveys

Secondary Dressing	Catalog / SKU#	HCPCS Code
Fabric Bandage (1 x 3)	A2002	N/A
Flex Tape (1 in.)	A60102	N/A
Knuckle Bandage (1 x 3.5)	A60103	N/A

AMERIGEL® HYDROGEL WOUND DRESSING

- Provides moist healing of dry wounds.
- Sustains moisture longer than standard hydrogels.
- Liquefies at body temperature for improved migration across all wound surfaces.
- Non-cytotoxic to healthy tissue.



*Podiatry Management Annual/Physicians Surveys

Podiatrists' #1 Choice for Wound/Ulcer Treatment and Treatment following Nail Surgery*

INDICATIONS:

For use as a wound dressing to manage pressure ulcers stages I-IV, stasis ulcers, diabetic skin ulcers, skin irritation, cuts, abrasions and skin irritations associated with peristomal care.

Reimbursable through Medicare and most private insurance plans.

Size	Catalog / SKU#	HCPCS Code
1 oz. Tube	A2001	A6248
3 oz. Tube	A20103	A6248

Saline Wound Wash (0.9% USP Normal Saline)

INDICATIONS:

To cleanse minor wounds, minor scrapes and minor burns.

Sterile

Drug, Preservative and Additive-FREE

Lockable Cap



AMERIGEL® SALINE WOUND WASH

- Painless sterile wound cleansing.
- Clinically appropriate pressure assists in removing debris and bacterial contaminants.
- No burning or stinging.
- Portable, 360° dispensing system sprays at any angle.

Size	Catalog / SKU#	HCPSC Code
4 oz.	A50114	N/A
7.1 oz.	A50117	N/A

AMERIGEL® HAND SANITIZER

- Three premium moisturizers to protect and hydrate your skin.
- Light fresh scent.
- Enjoy the confidence of using a trusted AMERIGEL® product.



Cleanse and Moisturize

INDICATIONS:

For handwashing when soap and water are not available.

Contains CDC-Recommended
70% alcohol

Size	Catalog / SKU#	HCPCS Code
6 oz. Tube	AS0106	N/A
16 oz. Pump Bottle	AS0116	N/A

100% Type 1 Bovine Collagen Effective in All Wound Phases

INDICATIONS:

Management of burns, sores,
blisters, ulcers
and other wounds.

Non-Hydrolyzed Tropocollagen

Biodegradable
and Biocompatible

No Additives, Synthetics or Fillers

Sterile

Rx Only

Reimbursable through Medicare
and most private insurance plans.



- Each gram absorbs **40-60x** its weight in fluid.
- Contributes up to **10x** more non-denatured collagen protein than competing products.
- Retained native helical structure enhances the vital “scaffolding” process for granulation tissue formation.

Size	Qty. per Box	Catalog / SKU#	HCPSC Code
1 gram	10	H40111	A6010



- Highly absorptive porous collagen sheet.
- Consists of fibrous triple-helix, tropocollagen structure.
- Can be used in combination with topical agents.
- Non-adherent.



**100% Type 1 Bovine Collagen
Effective in All Wound Phases**

INDICATIONS:

Management of burns, sores, blisters, ulcers and other wounds.

Non-Hydrolyzed Tropocollagen

Biodegradable
and Biocompatible

No Additives, Synthetics or Fillers

Sterile

Rx Only

Reimbursable through Medicare
and most private insurance plans.

Size	Qty. per Box	Catalog / SKU#	HCPCS Code
1" x 1"	10	H40220	A6021
2" x 2"	10	H40221	A6021
4" x 4"	10	H40224	A6021
7" x 7"	10	H40225	A6023

100% Type 1 Bovine Collagen Gel

INDICATIONS:

Management of abrasions, cuts, superficial injuries, severe sunburns, partial and full-thickness wounds, pressure ulcers, venous stasis ulcers, 1st and 2nd degree burns, ulcers caused by mixed etiologies, surgical wounds, donor sites and grafts.

Non-Hydrolyzed Tropocollagen
Biodegradable and Biocompatible
No Additives, Synthetics or Fillers
Sterile

Reimbursable through Medicare and most private insurance plans.



- Ideal for use when collagen powder or collagen sheets are not an option.
- Able to be used on any wound bed, including tunneling wounds.
- Can provide the benefits of collagen for treating 1st and 2nd degree burns.

Size	Qty. per Box	Catalog / SKU#	HCPSC Code
30 gram	1	H40401	A6011



CALCIUM ALGinate DRESSING

- Absorbs up to **20x** its weight in wound fluids.
- Quickly forms a hydrophilic gel to protect and cover the wound bed.
- Easy to pack, tuck or apply over irregularly shaped wounds.



Moderate to Heavy Draining Wounds

INDICATIONS:

Pressure ulcers, arterial ulcers, venous ulcers, diabetic ulcers, donor sites, post-operative wounds, dermal lesions, cuts & abrasions.

Sterile

No Natural Rubber Latex

Reimbursable through Medicare and most private insurance plans.



Size	Qty. per Box	Catalog / SKU#	HCPSC Code
2" x 2"	15	180121	A6196
4.25" x 4.25"	15	180125	A6197

Moderate to Heavy Draining Wounds

INDICATIONS:

Neuropathic ulcers, venous ulcers, open surgical wounds, donor sites, closed intact surgical wounds, lacerations, superficial burns, abrasions and skin tears.

Sterile

No Natural Rubber Latex

Reimbursable through Medicare and most private insurance plans.



- High performance polyurethane foam able to absorb large amounts of wound fluids.
- Breathable, self-adhesive border provides a waterproof barrier to outside contaminants.
- All-in-one dressing, ideal for use under compression.

Border Size	Pad Size	Qty. per Box	Catalog / SKU#	HCPSC Code
1" x 3.5"	1" x 1"	12	190220	A6413
4" x 4"	2" x 2"	12	190221	A6212
3" x 4"	2" x 2"	12	190222	A6212
6" x 6"	4" x 4"	12	190224	A6212

- Composed of breathable, non-adherent, high performance polyurethane foam.
- Able to absorb large amounts of wound fluids.
- Soft, conformable pads are ideal for use under compression.



Size	Qty. per Box	Catalog / SKU#	HCPSC Code
2" x 2"	12	190121	A6209
4.25" x 4.25"	12	190125	A6210

Moderate to Heavy Draining Wounds

INDICATIONS:

Neuropathic ulcers, venous ulcers, open surgical wounds, donor sites, closed intact surgical wounds, lacerations, superficial burns, abrasions and skin tears.

Sterile

No Natural Rubber Latex

Reimbursable through Medicare and most private insurance plans.

For Use as a Primary or Secondary Dressing

INDICATIONS:

Neuropathic ulcers, venous ulcers, open surgical wounds, donor sites, closed intact surgical wounds, lacerations, superficial burns, abrasions and skin tears.

Sterile

No Natural Rubber Latex

Reimbursable through Medicare and most private insurance plans.



AMER[®]
HEALTH CARE X

**BORDERED GAUZE
DRESSING**

- Flexible non-woven cloth border is hypoallergenic.
- Non-adherent cotton pads absorb wound exudate.
- Water-resistant to help protect the wound from contaminants.



Border Size	Pad Size	Qty. per Box	Catalog / SKU#	HCPSC Code
2" x 2"	1" x 1"	15	1G0220	A6219
4" x 4"	2" x 2"	15	1G0221	A6219
6" x 6"	4.25" x 4.25"	15	1G0226	A6220

AMER[®] **HEALTH CARE X** **RETENTION** **TAPE**

- Secure dressings and devices gently.
- Soft, comfortable tape uses gentle adhesive for sensitive skin.
- Breathable, non-woven dressing conforms to joints, bony prominences, and flexible body parts.
- Water resistant dressing provides dependable adherence in high-moisture prone areas.



For Use as a Secondary Dressing

RECOMMENDED FOR:

Securing primary dressings, catheters, feeding tubes, and similar devices.

No Natural Rubber Latex

Reimbursable through Medicare and most private insurance plans.

Size	Qty. per Box	Catalog / SKU#	HCPSC Code
2" x 11 Yards	1	1T0154	A4452



AMERX HEALTH CARE CORPORATION

164 Douglas Road, E
Oldsmar, FL 34677

AMERXHC.com
800.448.9599

ADVANTAGEN™ is a trademark of AMERX Health Care Corporation.
AMERIGEL®, AMERX®, HELIX3®, EXTREMIT-EASE®, Turn-Key DME® and Oakin® are
Registered Trademarks of AMERX Health Care Corporation.
Copyright© 2022 AMERX Health Care Corporation. All rights reserved.